

FillerCollector™

For gentle and effortless harvest of natural fat

In the past, fat harvesting was often a difficult, time-consuming process with unpredictable results.

But now, the all new FillerCollector™ simplifies the collection of injectable fat for autologous fat transfer. The FillerCollector™ is the most advanced fat harvesting system available for small to medium amounts of fat. The filterless, reusable collection device for fat harvesting, the FillerCollector™ is fast, easy and gentle on fat.

The FillerCollector™ is a better way to harvest fat.

Gentle

- Gentle on fat to help improve the chances of long-term fat survival*
- No centrifuging or other processing needed — saves time and protects fat cells from harm
- Pressure equalization inside container helps protect fat cells
- Gentle, passive blending of fat within container – improves fat consistency from multiple harvest sites

Effortless

- Simplified collection process
- Unique design continuously separates fat from rinsing fluid — no additional washing needed
- Easy fat extraction from the container
- Filterless — no clogging

Efficient

- Harvests high-quality fat**, which is ready for immediate transfer
- Pre-filter traps fibrous strains for more homogeneous fat
- Sterile, closed loop system
- Robust and durable components
- Autoclavable
- Harvests 100 ml in just 10 minutes

Gentle fat harvesting is the key to success.



human med

¹ Sasaki GH. Water-assisted liposuction for body contouring and lipoharvesting: safety and efficacy in 41 consecutive patients. *Aesthet Surg J.* 2011 Jan;31(1):76-88.

² K. Ueberreiter et al., BEAULI™ - A New and Easy Method for Large-Volume Fat Grafts. *Handchir Mikrochir Plast Chir* 2010; 42(6): 379-385.)

Ordering Information

Description		Units	REF
FillerCollector™		Complete set includes 1 FillerCollector™ (supplied non-sterile, autoclavable) 2 Disposable Sets FillerCollector™ (sterile silicone washer and plug), 2 Extraction Cannulas 300mm (sterile)	660000
Disposable Set FillerCollector™ (Each set consists of 2 silicone washers and 1 plug)	Sterile, for single use	5 sets	660023-5
Extraction Cannula 300mm	Sterile, for single use	5 sets	665010-5

Caution:

U.S. federal law restricts this device to sale by or on the order of a physician.

FillerCollector™ is a trademark owned by human med AG, Germany.

Patents pending.

©human med AG. 10/2013

REF 9001001

Contact Information

human med AG

Wilhelm-Hennemann-Str. 9 • 19061 Schwerin • Germany
Phone +49 (0) 385 / 395 70-0 • Fax +49 (0) 385 - 395 70-29
info@humanmed.com • www.humanmed.com

